For product information, contact:  
Lisätietoja tuotteesta antaa:  
För produktinformation, kontakta:  
Voor productinformatie, neem contact op met:  
Para información sobre el producto, contacte a:  
Pour obtenir des informations sur le produit, veuillez contacter :  
Produktinformationen:  
Per informazioni sui prodotti, contattate:

# Olsen MEDICAL®

3230 Commerce Center Place  
Louisville, KY 40211  
(502) 772-4280 <> (800) 297-6344

REF

Reference number / Tuotenumero / Referensnummer / Referentienummer / Número de referencia / Numéro de référence / Bezugsnummer / Numero di riferimento

LOT

Lot number / Eränumero / Satsnummer / Partijnummer / Número de lote / Numéro de lot / Seriennummer / Numero di lotto

See instructions for use / Lue käyttöohjeet / Se bruksanvisningen / Zie gebruiksaanwijzing / Veá las instrucciones de uso / Voir les instructions d'utilisation / Gebrauchsanweisung beachten / Consultare le istruzioni per le modalità d'uso

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No Latex

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

## ENGLISH

### DISPOSABLE ELECTROSURGICAL DEVICES

**Description:** Olsen Medical disposable electrosurgical devices are sold as sterile, single use devices, unless otherwise noted. These devices are designed to function as the active electrodes in electrosurgical systems. This family of devices includes, but is not limited to: electrosurgical pens, electrode handles, forceps, electrodes, cables, smoke evacuators, holster, tip cleaning pads, or various combination packages of these devices.

#### Handling:

- Inspect package for damage before use. Discard if damaged or seal is broken.
- Inspect device after opening in a sterile field. Discard if damaged.
- Device use and application to be determined by doctor.

Note: Remove and discard the tip protector on all electrodes before use.

**Storage:** Store at +5°C to +30°C. Relative humidity: 30% – 75%, non-condensing.

**Sterility:** Sterility is guaranteed up to the expiration date unless package is damaged or seal is broken.

**Sterilization:** Single Use Only. Dispose after use. Do not resterilize.

#### Caution:

- Damaged electrosurgical devices may cause electrical burns to the patient or doctor. Inspect every device before use.
- Active electrosurgical devices must be placed away from the patient when not in use. A holster is recommended.
- During use, keep the active electrode clean; build up of desiccated tissue will cause poor performance. A tip cleaning pad is recommended.
- Use with active electrodes only in RF surgical applications. Consult RF Generator owner's manual.
- Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the point of contact.
- Improperly seated electrode may cause burn or injury. Check that the electrode is compatible with the pencil or handle and ensure that the insulated sleeve fits securely inside the handle. Do not use electrode if it does not fit securely inside the handle.
- Do not bend, cut, or modify electrode before or during surgery. If bent, it may break the device or cause damage to insulation.
- Keep the voltage/power settings at the lowest possible level to achieve the desired effects.
- For Manually Extendable Electrocautery Apparatus, do not attempt to dislodge or remove extendable electrode to avoid breaking the device.

- After use, safely discard in accordance with established procedures for biohazardous waste.

**Warning:** Careless use of any element in electrosurgical systems may cause serious burn. The use of the device with safely functioning and compatible handles, cords, electrosurgical generators and accessories is the user's responsibility. Read and understand all precautions and instructions before attempting to use any electrosurgical device.

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