


BAYONET BIPOLAR TONSIL FORCEPS

Bayonet Tonsil Forceps, Bipolar, Insulated
Angled down, Smooth 1.5 mm Tips
7-1/2" (19.1cm), Working Length 3" (7.6cm)
Single Use, 10 per box



20-1381I Forceps

20-1381KI Kit including Forceps,
12' Cord (3.6M), and
Surgi-Wipe

Olsen MEDICAL[®]

Key Features

- Fully insulated with exposure only on the inside of the tips for increased safety
- Bipolar Coagulation for added security
- Offered as forceps alone or a convenient kit with cord and Surgi-Wipe
- Cost Effective
 - No processing costs
 - No repair or reinsulation
 - A fraction of the cost of a reusable
- Sterile and ready to use